**MODEL: UN300 A**

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Please read this manual thoroughly before use of the Device

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# CONTENTS

SUBJECTS	PAGE NO
1. GENERAL INFORMATION	5
2. SAFETY GUIDELINES	7
3. STARTING-UP/INSTALLATION	7
4. OPERATING PRINCIPLES	9
5. PREVENTIVE MAINTENANCE	11
6. TROUBLESHOOTING	11
7. TECHNICAL CHARACTERISTICS	11
8. LIST OF ACCESSORIES	11
9. GUARANTEE (TERMS & CONDITIONS)	12
10. CERTIFICATES	13
11. EMC DECLARATION	16
12. MANUFACTURER INFORMATION	18

## GLOSSARY OF SYMBOLS

	CAUTION CONSULT ACCOMPANYING DOCUMENTS		CLASS II EQUIPMENT		CONFORMITE EUROPE		NOT MOVED WITH HOOK
	DATE OF MANUFACTURE	<b>IPX1</b>	PROTECTION AGAINST VERTICALLY FALLING DROPS		THIS SIDE UP		SHOULD BE DISPOSED ACCORDING TO WEEE DIRECTIVES
	MANUFACTURE		TYPE B		FRAGILE, HANDLE WITH CARE		
<b>SN</b>	SERIAL NUMBER		TYPE BF		KEEP DRY		

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## 1.GENERAL INFORMATION

The device, NEBTIME ULTRASONIC NEBULIZER has been developed for aerosol & steam treatment of liquids.

In the conversion of liquids into vapour/mist (tiny particles of 0.5 – 6 µm), instead of heating, ultrasonic frequency technique has been applied. Frequencies produced by microprocessor based electronic system has been converted to mechanical vibrations with the help of Piezoceramic crystal which is responsible for the conversion of liquid placed into the flask of the device to vapour without heating. The evaporated liquid compel to come out of the flask reach to the patient with the help of debacterified environmental air pumped continuously into the flask by a fan fitted into the device. If require the out coming vapour can be warmed with the help of a heater placed at the outlet.

NEBTIME ultrasonic nebulizer has been manufactured to use the properties and characteristics

Microprocessor controlled

Plastic cabinet

Low electric consumption

Easy to replace Piezoceramic crystal

Drug application

Owned by it most effectively.

## FIELDS OF APPLICATION

Moisturizing of respiratory and narcotic gases,

Acute and chronic Bronchitis,

Bronchial Asthma,

Pneumonia, specific and nonspecific lung diseases,“

## SIDE EFFECTS

If the device is not cleaned regularly, bacteria which may form in the water reservoir, may affect the patient.

In case of using ventolin derivative drugs with cold steam, these drugs may have side effects.

## PATIENT POPULATION

Patients over 2 years of age and more than 10 kg with severe respiratory distress and treatment.

## CONTRAINDICATIONS

This product should not be used in patients with unstable and high blood pressure, cardiac irritability, increased heart rate, unconsciousness, or unable to breathe spontaneously.

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## 2.SAFETY GUIDELINES

Please go through the user guide A – Z before using the device.

We suggest to control the compatibility of electrical line voltage with the working voltage of the device along with electrical cord and plug to be used for the device. Before moving the device from one place to another please unplug the power supply cord.

Do not operate the device located very close to heat sources.

Do not put any item into the air filter vents of the device.

Keep the device away from intervention of children.

Please switch off the device using OFF/ON button and unplug the power supply cord when not in use.

Do not pour water on the device.

Consult with the troubleshooting section of this manual when you face any malfunctioning of the device.


Do not use any other liquids except from doctor's prescribed liquids in this device."

## 3. STARTING-UP / INSTALLATION

Please control the packing and the device on arrival. If there is any deformation or damage in the packing or in the inside materials, please inform the courier or the seller of the goods immediately. Take out from the package all components carefully controlling crash, breakage or any other damages.

**Let stand four (4) hours in user environment before start using the device.**

Place the device on plain and solid floor. Connect the power supply cord into the plug of AC 230V.

Fill the water flask of the device with cooled boiled or distilled water to the maximum level mark and close the lid and vapor tube connect to lid . Now the device is ready to use.

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## 4. OPERATING PRINCIPLES

Run the device operating OFF/ON switch to ON positing.

Set the vapour level as per requirement using NEB button of the device

Take special care keeping water level of the flask above minimum mark unless Piezo Crystal may damage. Water level is low , the device automatically stops and gives an audible alarm. Please , add enough water , before running on the device .

### **Device for use with medicine :**

put the medicine in to cup medicine cup attach on lid.  
close lid  
start device

Switch Off the device using OFF/ON switch when not in use and unplug the power supply cord.

**Keep the device away from children.**

Like most of electronic equipments NEBTIME too effects with Radio Frequency waves (Parasite). Please take care of not using mobile phone around the device.

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5 adet ilaç haznesi

## 5. PREVENTIVE MAINTENANCE

Please switch off the device & unplug power supply cord before maintenance work

### CLEANING OF OUTER CABINET

Do not use any liquid cleaning agent directly of the device. Using of petroleum based solvents and cleaning agents are prohibited. Clean the outer cabinet with pad wetted with soap water and well twisted.

### CLEANING OF WATER FLASK

Please clean the flask with well twisted pad wetted with soap water first and then wash the pad with plenty of water and wipe the flask several times.

## 6. TROUBLE SHOOTING

Please consult with authorized technical service if you face any unexpected situation with the device.

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## 9. GUARANTEE TERMS & CONDITIONS

- Guarantee period commences from the date of deliver of the product to client.
- In case of malfunctioning within the warranty period, repair time is added to the warranty period. Maximum repairing period is 30 (Thirty) working days.
  - If you have any objection regarding the device, please return the device to the company from where you bought the device or to the manufacturer by prepaid cargo along with the sales certificate and the guarantee certificate.
- Problems appeared due to wrong or misuse of the device is not covered under warranty.
- Opening of the device by unauthorized person will cause device to remain outside of the scope of the warranty.
- The guarantee for the device will only be valid once the consumer accepts after reading all terms and conditions mentioned here.
- Consumable parts are not covered under warranty

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# ELECTROMAGNETIC COMPATIBILITY (EMC) DECLARATION

This equipment generates, uses and can radiate radio frequency (RF) energy. This equipment may cause electromagnetic interference if it is not installed and used in accordance with the manual. This device has been tested for Medical Products in accordance with EN-60601-1-2:2015 and has been determined to comply with acceptable limits. These limits indicate that the device provides an acceptable level of protection against electromagnetic interference (EMC) if the device is used as specified in the manual. This device may be affected by portable and mobile RF communication devices. This equipment must not be stored together with other equipment. See the following tables for more information about this device and EMC.

## GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC EMISSIONS

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment

Emission test	Compliance	Electromagnetic environment – guidance
Emissions RF CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.  The device is suitable for use in all establishments including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Emissions RF CISPR 11	Class B	
Harmonic emissions used for domestic purposes. IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

## GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC EMISSIONS

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment

Immunity test	IEC 60601 Testlevel	Compliance level	Electromagnetic environment – guidance
Electrostatic discharges IEC 61000-4-2	± 8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	± 8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ±2 kV Input/output lines: ±1 kV 100 kHz repetition frequency	Power supply lines: ±2 kV Input/output lines: ±1 kV 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surges Line-to-line IEC 61000-4-5	± 0.5 kV, ± 1 kV	± 0.5 kV, ± 1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surges Line-to-ground IEC 61000-4-5	± 0.5 kV, ± 1 kV, ± 2 kV	± 0.5 kV, ± 1 kV, ± 2 kV	
Voltage dips IEC 61000-4-11	0% UT; 0.5 cycle 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°  0% UT; 1 cycle and 70% UT; 25/30 cycle Single phase: 0°	0% UT; 0.5 cycle 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°  0% UT; 1 cycle and 70% UT; 25/30 cycle Single phase: 0°	Mains power quality should be that of a typical commercial or hospital environment.
Voltage interruptions IEC 61000-4-11	0% UT; 250/300 cycle	0% UT; 250/300 cycle	
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Not: UT is the a.c. mains voltage prior to application of the test level.

Conducted RF IEC 61000-4-6	150kHz to 80MHz: 3Vrms 6 Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	150kHz to 80MHz: 3Vrms 6 Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	Portable and mobile RF communications equipments should be used no closer to any part of the device including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 0.35\sqrt{P}$ $d = 1.2\sqrt{80 \text{ MHz ile } 800 \text{ MHz}}$ $d = 1.2\sqrt{800 \text{ MHz ile } 2.7 \text{ GHz}}$ $d = 2.3\sqrt{P}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	10V/m 80% Am at 1kHz	10V/m 80% Am at 1kHz	

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Note 1: 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

c. Calibration for current injection clamps shall be performed in a 150 Ω system.

d. If the frequency stepping skips over an ISM or amateur band, as applicable, an additional test frequency shall be used in the ISM or amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.

e. Applicable to ME EQUIPMENT and ME SYSTEMS with RATED input current less than or equal to 16 A / phase and ME EQUIPMENT and ME SYSTEMS with RATED input current greater than 16 A / phase.

Recommended separation distances between portable and mobile RF communications equipment and the device

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter(m)			
	150 kHz to 80 MHz (out ISM and amateur radio bands) d= 1.2 √ P	150 kHz to 80 MHz (out ISM and amateur radio bands) d= 0.6 √ P	80 MHz to 800 MHz d= 1.2 √ P	800 MHz to 2.7 GHz d= 2.3 √ P
0,01	0,12	0,06	0,12	0,23
0,1	0,38	0,19	0,38	0,73
1	1,2	0,6	1,2	2,3
10	3,8	1,9	3,8	7,3
100	12	6	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.  
Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

Note 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Radiated RF	Test Frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Modulation b) (W)	Distance (m)	Immunity Test Level(V/m)
IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipments)	385	380-390	TETRA 400	Pulse modulation b) 18 Hz	1,8	0,3	27
	450	380-390	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0,3	28
	710	704-787	LTE BAND1 13, 17	Pulse modulation b) 217 Hz	0,2	0,3	9
	745						
	780						
	810	800-960	GSM 800/900, TETRA 800, IDEN 820, LTE Bant 5	Pulse modulation b) 18 Hz	2	0,3	28
	870						
930							
1720	1700-1990	GSM 1800, CDMA 1900, GSM 1900,	Pulse modulation b) 217 Hz	2	0,3	28	
1845		DECT, LTE Bant1, 3,4,25, UMTS					
1970		Bluetooth, WLAN 802.11b/g/n, RFID 2450, LTE Bant 7	Pulse modulation b) 217 Hz	2	0,3	28	
2450	2400-2570						
5240	5100-5800	WLAN 802.11a/n	Pulse modulation b) 217 Hz	0,2	0,3	9	
5240							
5785							

Note: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

a. For some services, only the uplink frequencies are included.

b. The carrier shall be modulated using a 50 % duty cycle square wave signal.

c. As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation:

$$E = 6 / d \sqrt{P}$$

Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.

## 12. ÜRETİCİ FİRMA / MANUFACTURER

### ELMASLAR

İmalat Tibbi Cihazlar İthalat İhracat A.Ş.

ASO 1. Org. San. Bölgesi Babürşah Cad. No: 17

Sincan /ANKARA/ TURKIYE

Tel: +90 312 394 80 01 Fax: +90 312 394 80 04

info@elmaslarmedikal.com.tr • elmaslarmedikal@gmail.com

www.elmaslarmedikal.com.tr

## **ELMASLAR A.S.**

ASO 1.OSB BABURSAH CD. NO: 17 SINCAN / ANKARA / TURKEY

Tel: +90 312 394 80 01 Fax: +90 312 394 80 04

info@elmaslarmedikal.com.tr • elmaslarmedikal@gmail.com

[www.elmaslarmedikal.com.tr](http://www.elmaslarmedikal.com.tr)